




PRODUCT SELECTION FOR MEDICAL PROVIDERS



PowerStep®



Unsupported Foot
tends to roll inward. The arch flattens, stretching the tissue while bones misalign.



Foot with Powerstep®
provides foot posture and support to correct bone alignment and prevent pain.

HOW DO POWERSTEP® ORTHOTICS WORK?

Powerstep® orthotics were developed by a podiatrist to provide a unique foot support system that goes beyond “foot comfort.” Powerstep® orthotics support the foot with a polypropylene, semi-flexible arch support that is encapsulated between two layers of advanced foam cushioning. An open arch area fully supports the foot and arch with enough spring and flexibility to maintain proper foot function and alignment.

Powerstep® products are available in a variety of styles, each with different levels of cushioning and support to meet the unique needs of your patients. Products are available for various levels of athletic activity, shoe styles and personal preferences.

POWERSTEP® ORTHOTICS ARE *PROVEN EFFECTIVE* IN CLINICAL STUDIES.

- Powerstep® is clinically proven to improve foot function and decrease foot pain.²
- As compared to walking with no orthosis, Powerstep® resulted in a 33.19% reduction in Plantar Fascia strain.⁷
- When compared to two other over-the-counter devices on the market, only Powerstep® exhibited significant decreases in peak rear-foot eversion.⁷
- Powerstep® orthotics are an effective, efficient and economical alternative to custom orthotics in treating plantar heel pain.³
- Orthotics incorporating a rigid, plastic component are superior in reducing pain associated with Plantar Fasciitis.⁴
- Foot orthotics may help prevent and treat lower back pain.⁵
- Foot orthotics are effective in treating ankle instability.⁶

ADVANTAGES OF IN-OFFICE DISPENSING OF PREFABRICATED ORTHOTICS:

- ✓ Affordable
- ✓ First-line therapy for mild conditions
- ✓ Immediate relief while waiting for custom orthotics
- ✓ An extra device for a second pair of shoes
- ✓ Training for custom orthotics or testing the need for them

Proven Clinically Effective*



Why Offer Powerstep® therapeutic products to your patients?

- **Powerstep® is the #1 Podiatrist-Recommended Brand of Prefabricated Orthotics.**¹
- **Powerstep® is clinically proven to improve foot function and decrease foot pain.**²
- **Increase Patient Satisfaction** with a true podiatrist-designed product that relieves heel and arch pain, is widely endorsed and referred by physicians and is accepted by the American Podiatric Medical Association.
- **In-Office Dispensing** allows you to offer patients immediate and affordable pain relief and first-line therapy for mild conditions while waiting for custom orthotics, if necessary.
- **Simple Solutions for Common Conditions**, easy-to-use-and-fit prefabricated products for effective pain relief of many common foot conditions.



PREFABRICATED ORTHOTICS

Exclusive to Medical Professionals*

- 2 | ProTech Control Maximum Support**
Control Full Length, Control Met, Control Wide, Control 3/4
- 3 | ProTech Full Support**
ProTech Full Length, Classic Plus, ProTech Met, ProTech 3/4
- 3 | ProTech CustomPost®**
- 3 | Diabetic Insoles**
Bi-Laminate, Tri-Laminate
- 4 | Retail Insoles**
Pinnacle Maxx, Wide Fit, Pinnacle, Pinnacle Pink, Pinnacle Plus, Original
- 4 | Cushioning Insoles**
Memory Foam, ComfortLast®
- 5 | Low Profile Insoles**
SlimTech®, SlenderFit®, Signature 3/4
- 5 | Toddler and Youth Styles**
PowerKids®, Pinnacle Junior, KidSport®

- 6 | CUSTOM ORTHOTICS**
Light Orthotics® by Powerstep

- 7 | SHOE ACCESSORIES**
Orthotic Modifiers
Orthotics Modification Kit

- 8-9 | SUPPORTIVE FOOTWEAR**
Recovery Sandals & Slides
Orthotic Slippers

- 10-11 | SUPPORTS & BRACES**
Plantar Fasciitis Supports, Sport Supports, Post-Op Shoe, Dorsal Night Splint

- 11 | COMPRESSION**
Performance Calf Sleeves, Compression Socks

- 12-13 | Merchandising Support**

*PROTECH DISPENSING

We have developed the Powerstep® ProTech Line as exclusive products distributed ONLY by medical professionals and for use ONLY by their patients. Accordingly, ProTech products cannot be displayed or sold through websites, catalogs, stores or any other sales arrangement that is accessible to the general public. Powerstep® products purchased from unauthorized resellers are not covered by our total satisfaction guarantee.

¹ Podiatry Management Online. "Which Brand of Pre-Fabricated Orthotics Do You Dispense in Your Office?" Survey. <http://podiatrym.com/photo.cfm?photoid=8884>.

² Springett, Otter and Barry. "A Clinical Longitudinal Evaluation of Pre-Fabricated, Semi-Rigid Foot Orthoses Prescribed to Improve Foot Function." *The Foot*. 17 (2007) 14-189.

³ Ring and Otter. "Clinical Efficacy and Cost Effectiveness of Bespoke Compared with Over-the-Counter Foot Orthoses for Plantar Heel Pain." *Journal of Foot and Ankle Research*. 2010 3 (Suppl 1): 022.

⁴ Hume, P, Hopkins W, Rome K, Maulder P, Coyle G, and Nigg B. "Effectiveness of Foot Orthoses for Treatment and Prevention of Lower Limb Injuries." *Sports Medicine*. 2008;38(9):759-79.

⁵ Menz HB, Dufour AB, Riskowski JL, Hillstrom HJ and Hannan MT. "Foot Posture, Foot Function and Low Back Pain: The Framingham Foot Study." *Rheumatology* (Oxford). Sept 17, 2013.

⁶ Douglas R Jr. "Chronic Ankle Instability: Can Orthotics Help?" *Podiatry Today*. October 2006. Vol 19-Iss 10.

⁷ Ferber, Reed and Hettinga, BA. "A Comparison of Different Over-the-Counter Foot Orthotic Devices on Multi-Segment Foot Biomechanics." *Prosthetics and Orthotics International*. May 26, 2015.

PROTECH EXCLUSIVE STYLES

Exclusive to Medical Professionals

ProTech	Main Feature	Pronation	Cushioning	Neuromas	Recommended Shoe Style
ProTech Full Length	Full length cushioning and support	Mild	Maximum		Athletic, Casual
ProTech Met	Support with built-in metatarsal pad to relieve ball-of-foot pain	Mild	Maximum	✓	Athletic, Dress Casual
ProTech Classic Plus	Full length support in a thinner profile	Mild	Moderate		Dress, Casual, Low profile shoes
ProTech 3/4	Slimmer 3/4 length support ideal for tighter fitting shoes	Mild	Moderate		Athletic, Dress, Casual Tight-fitting shoes, Cleats
ProTech CONTROL	Main Feature	Pronation	Cushioning	Neuromas	Recommended Shoe Style
ProTech Control	Full length cushioning and support	Moderate to Severe	Maximum		Athletic, Casual
ProTech Control Met	Support with built-in metatarsal pad to relieve ball-of-foot pain	Moderate to Severe	Maximum	✓	Athletic, Casual
ProTech Control Wide	Wider toe box and greater support to prevent bottoming-out under weight	Moderate to Severe	Maximum		Athletic, Casual, Wide-width and Therapeutic footwear
ProTech Control 3/4	Slimmer 3/4 length support ideal for tighter fitting shoes	Moderate to Severe	Moderate		Athletic, Dress, Casual Tighter-fitting shoes, Cleats



PROTECH CONTROL



PROTECH CONTROL MET
with built-in metatarsal pad



PROTECH CONTROL WIDE
with wider profile



PROTECH CONTROL 3/4

1015-01	PROTECH CONTROL	Men's 3-16+ Women's 5-12	Maximum, dual-layer cushioning from heel to toe with firmer arch and 2° heel post for proper alignment and greater correction.
1019-01	PROTECH CONTROL MET	Men's 4-15 Women's 6-12	Firmer arch support with plush dual-layer cushioning and built-in metatarsal pad to spread and cushion the met heads. 2° heel post for improved foot alignment. Provides relief for Morton's neuroma and metatarsalgia.
1017-01	PROTECH CONTROL WIDE	Men's 4-15 Women's 6-12	Firmer arch and 2° heel post in a wider profile for those needing greater support and correction or have a wider foot. Ideal for wide-width footwear. Fits 3E-6E.
1015-03	PROTECH CONTROL 3/4	Men's 4-13½ Women's 6-15½	Firmer arch support and 2° heel post with moderate cushioning in a 3/4 length leaves extra room at the toes. Fits a variety of dress, casual and sports shoes.

Sizing	Asm	A	B	C	D	E	F	G	H	I/J	K	L
Men's Shoe	3 – 3½	4 – 4½	5 – 5½	6 – 6½	7 – 7½	8 – 8½	9 – 9½	10 – 10½	11 – 11½	12 – 13	14 – 15	16+
Women's Shoe	5 – 5½	6 – 6½	7 – 7½	8 – 8½	9 – 9½	10 – 10½	11 – 11½	12	—	—	—	—
CONTROL FULL	1015-01AS	1015-01A	1015-01B	1015-01C	1015-01D	1015-01E	1015-01F	1015-01G	1015-01H	1015-01I	1015-01K	1015-01L
CONTROL MET	—	1019-01A	1019-01B	1019-01C	1019-01D	1019-01E	1019-01F	1019-01G	1019-01H	1019-01I	1019-01K	—
CONTROL WIDE	—	1017-01A	1017-01B	1017-01C	1017-01D	1017-01E	1017-01F	1017-01G	1017-01H	1017-01I	1017-01K	—

Sizing	Asm	A	B	C	D	E	F
Men's Shoe	4 – 4½	5 – 6½	7 – 8½	9 – 9½	10 – 10½	11 – 11½	12 – 13½
Women's Shoe	6 – 6½	7 – 8½	9 – 10½	11 – 11½	12 – 12½	13 – 13½	14 – 15½
CONTROL 3/4	1015-03AS	1015-03A	1015-03B	1015-03C	1015-03D	1015-03E	1015-03F

ProTech Control styles are covered under US Patent D661,076 and are HCPCS L3040 Approved.



PROTECH EXCLUSIVE STYLES, DIABETIC INSOLES



PROTECH FULL LENGTH



PROTECH MET
with built-in metatarsal pad



PROTECH CLASSIC PLUS
for a thinner profile



PROTECH 3/4

1005-01	PROTECH FULL LENGTH	Men's 3-16+ Women's 5-12	Maximum, dual layer cushioning from heel to toe with semi-rigid arch support for proper alignment and improved foot function. Total contact support. US Design Patent D846,853
1009-01	PROTECH MET	Men's 4-15 Women's 6-12	Semi-rigid arch support, plush dual-layer cushioning and built-in met pad to spread and cushion the met heads. Relief for Morton's neuroma and metatarsalgia. US Design Patent D846,853
1003-01	PROTECH CLASSIC PLUS	Men's 3-15 Women's 5-12	Moderate, dual-layer cushioning with semi-rigid arch support in a thinner profile to better fit a variety of shoes. US Design Patent D846,853
1005-03	PROTECH 3/4	Men's 4-13½ Women's 6-15½	Arch support and moderate cushioning in a 3/4 length leaves extra room at the toes. Fits a variety of dress, casual and sports shoes.

Sizing	Asm	A	B	C	D	E	F	G	H	I/J	K	L
Men's Shoe	3 - 3½	4 - 4½	5 - 5½	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13	14 - 15	16+
Women's Shoe	5 - 5½	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12	—	—	—	—
PROTECH FULL	1005-01AS	1005-01A	1005-01B	1005-01C	1005-01D	1005-01E	1005-01F	1005-01G	1005-01H	1005-01I	1005-01K	1005-01L
CLASSIC PLUS	1003-01AS	1003-01A	1003-01B	1003-01C	1003-01D	1003-01E	1003-01F	1003-01G	1003-01H	1003-01I	1003-01K	—
PROTECH MET	—	1009-01A	1009-01B	1009-01C	1009-01D	1009-01E	1009-01F	1009-01G	1009-01H	1009-01I	—	—

Sizing	Asm	A	B	C	D	E	F
Men's Shoe	4 - 4½	5 - 6½	7 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13½
Women's Shoe	6 - 6½	7 - 8½	9 - 10½	11 - 11½	12 - 12½	13 - 13½	14 - 15½
ProTech 3/4	1005-03AS	1005-03A	1005-03B	1005-03C	1005-03D	1005-03E	1005-03F

ProTech styles are HCPCS L3040 approved.



PROTECH CUSTOMPOST®



1020-01	PROTECH CUSTOMPOST®	Men's 4-14½ Women's 6-12
---------	----------------------------	-----------------------------

Includes two sets of snap-on heels for heel posting and heel heights. Snap-on Neutral, 2° and 4° heel post options can be applied for either varus or valgus posting. Snap-on heel heights, 1/8" and 1/4". Firmer, medical grade support shell. Dual-layer, shock absorbing EVA foam cushioning for comfort. US Patent 7,913,429

Sizing	Asm	A	B	C	D	E	F	G
Men's Shoe	4 - 4½	5 - 6½	7 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13½	14 - 14½
Women's Shoe	6 - 6½	7 - 8½	9 - 10½	11 - 11½	12	—	—	—
CUSTOMPOST	1020-01AS	1020-01A	1020-01B	1020-01C	1020-01D	1020-01E	1020-01F	1020-01G



BI-LAMINATE

TRI-LAMINATE

2002-01	BI-LAMINATE DIABETIC INSOLES	3 PAIR/ PACK	Men's 6-15 Women's 5-12
2003-01	TRI-LAMINATE DIABETIC INSOLES	3 PAIR/ PACK	Men's 6-15 Women's 5-12

Heat formable 20 Shore A EVA top cover layer for total contact comfort and support and firm 40 Shore A EVA base. Tri-Laminate insole includes PORON® Cushioning layer for shock absorption. Lightweight, durable and heat formable.

HCPCS A5512 Approved

Sizing	X-Small	Small	Medium	Large	X-Large	XX-Large
Men's Shoe	—	6 - 7	8 - 9	10 - 11	12 - 13	14 - 15
Women's Shoe	5 - 7	8 - 9	10 - 11	12	—	—
BI-LAMINATE	2002-01XS	2002-01S	2002-01M	2002-01L	2002-01XL	2002-01XXL
TRI-LAMINATE	2003-01XS	2003-01S	2003-01M	2003-01L	2003-01XL	2003-01XXL

RETAIL INSOLES, CUSHIONING INSOLES



PINNACLE MAXX
firmer with heel posting



WIDE FIT
firmer with wider toe box



PINNACLE



PINNACLE PINK

5015-01	PINNACLE MAXX	Men's 3-16+ Women's 5-12	Firmer arch with a slightly-angled, 2° heel post for greater correction of moderate overpronation. Dual-layered cushioning. US Patent D661,076
5017-01	WIDE FIT	Men's 3-16+ Women's 5-12	Firmer arch and dual-layer cushioning. Wider toe box fits widths 3E-6E. Suitable for those needing heavier-weight support or a wider foot.
5005-01	PINNACLE	Men's 3-16+ Women's 5-12	Dual-layer cushioning with full support polypropylene arch shell and neutral heel. US Design Patent D846,853
5006-01	PINNACLE PINK	Women's 6-11½	Dual-layer cushioning with full support polypropylene arch shell and neutral heel. US Design Patent D846,853



PINNACLE PLUS
with built-in met pad



ORIGINAL
with thinner cushioning



MEMORY FOAM
with memory foam cushion



COMFORTLAST®
with formed arch and gel heel

5009-01	PINNACLE PLUS	Men's 4-15 Women's 6-12	Dual-layer cushioning with built-in met pad to position the met heads to relieve pain from neuromas. Full support polypropylene arch shell and neutral heel. US Design Patent D846,853
5001-01	ORIGINAL	Men's 3-16+ Women's 5-12	Moderate cushioning and full arch support in a thin, low-profile design to fit more shoe types. Neutral heel. US Design Patent D846,853
5018-01	MEMORY FOAM	Men's 4-15 Women's 6-12	Plush, slow-recovery memory foam with full arch support. Ideal for sensitive feet, fat pad atrophy and pressure relief. US Design Patent D846,853
5002-01	COMFORTLAST®	Men's 3-14½ Women's 5-12	Foam throughout with shock-absorbing gel at heel. Contoured shape with anatomical gel arch support built-in. No polypropylene shell.

Sizing	Asm	A	B	C	D	E	F	G	H	I/J	K	L
Men's Shoe	3 - 3½	4 - 4½	5 - 5½	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13	14 - 15	16+
Women's Shoe	5 - 5½	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12	—	—	—	—
PINNACLE MAXX	5015-01AS	5015-01A	5015-01B	5015-01C	5015-01D	5015-01E	5015-01F	5015-01G	5015-01H	5015-01I	5015-01K	5015-01L
WIDE FIT	5017-01AS	5017-01A	5017-01B	5017-01C	5017-01D	5017-01E	5017-01F	5017-01G	5017-01H	5017-01I	5017-01K	5017-01L
PINNACLE	5005-01AS	5005-01A	5005-01B	5005-01C	5005-01D	5005-01E	5005-01F	5005-01G	5005-01H	5005-01I	5005-01K	5005-01L
PINNACLE PINK	—	5006-01A	5006-01B	5006-01C	5006-01D	5006-01E	5006-01F	—	—	—	—	—
PINNACLE PLUS	—	5009-01A	5009-01B	5009-01C	5009-01D	5009-01E	5009-01F	5009-01G	5009-01H	5009-01I	5009-01K	—
ORIGINAL	5001-01AS	5001-01A	5001-01B	5001-01C	5001-01D	5001-01E	5001-01F	5001-01G	5001-01H	5001-01I	5001-01K	5001-01L
MEMORY FOAM	—	5018-01A	5018-01B	5018-01C	5018-01D	5018-01E	5018-01F	5018-01G	5018-01H	5018-01I	5018-01K	—

Sizing	1	2	3	4	5	6
Men's Shoe	3 - 4½	5 - 6½	7 - 8½	9 - 10½	11 - 12½	13 - 14½
Women's Shoe	5 - 6½	7 - 8½	9 - 10½	11 - 12	—	—
COMFORTLAST®	5002-011	5002-012	5002-013	5002-014	5002-015	5002-016



5020-03	SIGNATURE DRESS 3/4	Men's 4-13½ Women's 6-15½	Semi-rigid arch support with soft, supple top cover fabric for stylish wear. 3/4 length ideal for shoes with less room at the toes. US Patents
5012-03	SLIMTECH® 3/4	Men's 4-13½ Women's 6-15½	Ideal for tighter-fitting shoes or cleats with less room at the toes. Exposed arch shell with non-slip pad at heel to hold in place. US Patent D860,617
5010-03	SLENDERFIT®, BLACK	Women's 5-10½	Ultra-slim tapered design for support in women's fashion footwear. Cushioning under ball of foot for increased comfort and protection. US Patent D718,025 & D748,387
5011-03	SLENDERFIT®, NATURAL	Women's 5-10½	New Natural color for more discreet wear. Ultra-slim tapered design for support in women's fashion footwear. Cushioning under ball of foot for increased comfort and protection. US Patent D718,025 & D748,387



5040-03	POWERKIDS® 3/4	Toddler 8½ - Youth 4	Firm contoured EVA arch support and deep heel cradle with medial 2-degree heel post for motion control. No polypropylene shell. US Patent D701,028
5041-01	PINNACLE JUNIOR	Toddler 11½ - Youth 6	Our popular Pinnacle support now available in Toddler and Youth sizes. Full length cushioning and semi-rigid polypropylene arch support shell. US Design Patent D846,853
5043-01	KIDSPORT®	Toddler 9 - Youth 6	Dual-layer cushioning from heel to toe. No polypropylene shell. Ideal for children needing improved comfort in shoes. US Design Patent D846,853

Sizing	As	A	B	C	D	E	F
Men's Shoe	4 - 4½	5 - 6½	7 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13½
Women's Shoe	6 - 6½	7 - 8½	9 - 10½	11 - 11½	12 - 12½	13 - 13½	14 - 15½
SIGNATURE 3/4	5020-03AS	5020-03A	5020-03B	5020-03C	5020-03D	5020-03E	5020-03F
SLIMTECH®	5012-03AS	5012-03A	5012-03B	5012-03C	5012-03D	5012-03E	5012-03F

Sizing	XS	SM	MD	LG
Women's Shoe	5 - 6	6½ - 7½	8 - 9	9½ - 10½
SLENDERFIT® Black	5010-03XS	5010-03S	5010-03M	5010-03L
SLENDERFIT® Natural	5011-03XS	5011-03S	5011-03M	5011-03L

Sizing	KB	KC	KD	KE	KF	KG
Toddler/Youth	T 8½ - 9½	T 10 - 11	T 11½ - 12½	T 13 - Y1	T 1½ - 2½	T 3 - 4
POWERKIDS®	5040-03KB	5040-03KC	5040-03KD	5040-03KE	5040-03KF	5040-03KG

Toddler/Youth	T9	T10	T11	T11 ½	T12	T13	Y1	Y2	Y3	Y4	Y5	Y6
PINNACLE JR	—	—	—	5041-01115	5041-0112	5041-0113	5041-011	5041-012	5041-013	5041-014	5041-015	5041-016
KIDSPORT®	5043-019	5043-0110	5043-0111	5043-0112	5043-0112	5043-0113	5043-011	5043-012	5043-013	5043-014	5043-015	5043-016

5 MINUTES TO RELIEF

Shapes a flexible material into a custom orthotic using light-curing technology

Patient leaves office with custom orthotics in same visit



Learn more at <https://www.lightorthotics.com>



Size	A	B	C	D	E	F	G	H	I/J	K	L
Fits Men	4 - 4½	5 - 5½	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12 - 13	14 - 15	16+
Fits Women	6 - 6½	7 - 7½	8 - 8½	9 - 9½	10 - 10½	11 - 11½	12	—	—	—	—
PRETHOTICS	1075-01A	1075-01B	1075-01C	1075-01D	1075-01E	1075-01F	1075-01G	1075-01H	1075-01I	1075-01K	1075-01L

Also available in 3/4 length and toddler/youth sizes












To learn more about Powerstep® Light Orthotics, contact your Sales Rep at 888-237-3668.

SHOE ACCESSORIES & MODIFICATIONS

Powerstep® shoe accessories can be used in conjunction with Powerstep® orthotics and Light Orthotics® for enhanced support.

*Bulk pack items are not packaged for individual resale.

SHOE ACCESSORIES & ORTHOTIC MODIFIERS

ARCH BOOSTERS		5090-01, 5090-06 Sold as Pair or in Bulk pack*	SM, REG	Can be worn with Powerstep® insoles to boost the support or applied directly to shoe for arch support. Ideal for severely flat feet/fallen arches. Adhesive backing.
BUNION FELT PADS		5099-03 Sold in Bulk pack*	UNIV	Designed to offset bunion pressure and cushion and protect the area from rubbing and chaffing. Adhesive backing.
ADJUSTABLE HEEL LIFT		5093-01 Sold as Each	SM-LG	Provides 3/8" additional heel height or peel away layers for 1/4" or 1/8" additional height. Polyester top cover and tapered edge. Ideal for leg length discrepancy.
CORK HEEL LIFT		5094-01 Sold as Each	Sizes 1-4	Composite cork construction with polyester top cover. Tapered edge. Available in 4mm, 6mm and 8mm heights. Ideal for leg length discrepancy.
HEEL LIFTS		5092-03 1/4" 5092-13 1/8" Sold in Bulk pack*	SM, LG	Provides 1/4" or 1/8" additional heel height to correct leg length discrepancy. Dense, durable EVA with adhesive backing.
FOREFOOT WEDGES		5091-03 Right 5091-13 Left Sold in Bulk pack*	SM, LG	Designed for placement under first ray for forefoot varus/forefoot supinatus. Made of durable EVA. Adhesive backing.
2° HEEL WEDGES		5095-01, 5095-06 Sold as Pair or in Bulk pack*	SM, REG	Provides two degrees of additional correction to Powerstep insoles when applied to the heel platform. Dense, durable EVA with adhesive backing.
4° HEEL WEDGES		5095-16 Sold in Bulk pack*	SM, REG	Provides four degrees of additional correction to Powerstep insoles when applied to the heel platform. Dense, durable EVA with adhesive backing.
IPK BALL OF FOOT CUSHIONS		5096-01, 5096-06 Sold as Pair or in Bulk pack*	UNIV	Hollowed center to offload pressure from dropped metatarsal head. Cushions and surrounds Intractable Plantar Keratoma (IPK) or deep calluses. Adhesive backing.
METATARSAL PADS		5097-01, 5097-06 Sold as Pair or in Bulk pack*	SM, REG	Provide extra cushioning to spread the metatarsal bones to relieve forefoot pain. Left and right design for a more exact fit. Adhesive backing.
METATARSAL TEARDROP PADS		5098-06 Sold in Bulk pack*	SM, REG	Help spread the metatarsal bones to alleviate pain attributed to metatarsalgia and Morton's neuroma. Adhesive backing.

5000KIT ORTHOTICS MODIFICATION KIT

Contains a collection of accessory items for practitioners to modify and enhance Powerstep® prefabricated orthotics. The included assortment, containing 336 individual pieces, includes the most common modifiers used to customize orthotics for common foot conditions, all packaged neatly in a sturdy, organized carrying case. Each item is available in convenient, affordable reorder bulk packs for restocking.

The Kit Includes:

- **24 pair Arch Boosters (SM, REG)** for extra support under arch and midfoot
- **48 Forefoot Wedges (SM, REG, Right, Left)** for forefoot varus/forefoot supinatus
- **12 pair 1/4" Heel Lifts & 12 pair 1/8" Heel Lifts** for leg length discrepancy
- **24 pair 2° Heel Wedges & 24 pair 4° Heel Wedges** to correct pronation/foot alignment
- **6 pair IPK Ball of Foot Cushions** for intractable plantar keratosis
- **12 Bunion Felt Pads** to cushion and offset pressure on bunions
- **12 pair Metatarsal Cushions (SM, REG) & 24 pair Metatarsal Teardrop pads (SM, REG)** to spread and cushion metatarsal bones





MEN'S FUSION SLIPPER
FUSION | 8000-30
Brown/Cream

New!
MEN'S FUSION SLIPPERS



MEN'S FUSION SLIPPER
FUSION | 8000-20
Charcoal/Gray

- Built-in arch and heel
- cup for support and
- stability. Faux fur lining
- and suede-like outer with
- textile bottom. Available
- in Charcoal or Brown.



MEN'S RECOVERY SANDAL
FUSION | 8000
Black/Gray



MEN'S RECOVERY SANDAL
FUSION | 8000
Brown/Tan



MEN'S RECOVERY SLIDE
FUSION | 8100
BLACK

WOMEN'S SUPPORTIVE FOOTWEAR



WOMEN'S RECOVERY SANDAL
FUSION | 8500-10
Black/Gray




WOMEN'S RECOVERY SANDAL
FUSION | 8500-40
Charcoal/Plum/Gray



WOMEN'S RECOVERY SANDAL
FUSION | 8500-35
Khaki/Coral/Tan



WOMEN'S RECOVERY SANDAL
FUSION | 8500-30
Brown/Coral/Tan



MENS SIZING	7	8	9	10	11	12	13	14
RECOVERY SANDAL BLACK	8000-1007	8000-1008	8000-1009	8000-1010	8000-1011	8000-1012	8000-1013	8000-1014
RECOVERY SANDAL BROWN	8000-3007	8000-3008	8000-3009	8000-3010	8000-3011	8000-3012	8000-3013	8000-3014
RECOVERY SLIDE	—	8100-1008	8100-1009	8100-1010	8100-1011	8100-1012	8100-1013	8100-1013
FUSION SLIPPER BROWN/CREAM	—	8800-3008	8800-3009	8800-3010	8800-3011	8800-3012	8800-3013	8800-3014
FUSION SLIPPER CHARCOAL/GRAY	—	8800-2008	8800-2009	8800-2010	8800-2011	8800-2012	8800-2013	8800-2014



WOMEN'S SIZING	6	7	8	9	10	11
BLACK/LIGHT GRAY	8500-1006	8500-1007	8500-1008	8500-1009	8500-1010	8500-1011
CHARCOAL/PLUM/GRAY	8500-4006	8500-4007	8500-4008	8500-4009	8500-4010	8500-4011
KHAKI/CORAL/TAN	8500-3506	8500-3507	8500-3508	8500-3509	8500-3510	8500-3511
BROWN/CORAL/TAN	8500-3006	8500-3007	8500-3008	8500-3009	8500-3010	8500-3011
CHARCOAL/TEAL/GRAY	8500-5006	8500-5007	8500-5008	8500-5009	8500-5010	8500-5011
RECOVERY SLIDE	8150-1006	8150-1007	8150-1008	8150-1009	8150-1010	8150-1011
LUXE® SLIPPER BLACK/CHARCOAL	8850-1006	8850-1007	8850-1008	8850-1009	8850-1010	8850-1011
LUXE® SLIPPER TAUPE/IVORY	8850-3006	8850-3007	8850-3008	8850-3009	8850-3010	8850-3011



WOMEN'S RECOVERY SANDAL
FUSION | 8500-50
Charcoal/Teal/Gray



WOMEN'S RECOVERY SLIDE
FUSION | 8150-10
Black



WOMEN'S LUXE® SLIPPERS
LUXE® | 8850-10
Black/Charcoal



WOMEN'S LUXE® SLIPPERS
LUXE® | 8850-30
Taupe/Ivory

PLANTAR FASCIITIS RELIEF, SUPPORTS & BRACING

SQUARE TOE POST-OP SHOE

7500-01 Each, SM-XL

Provides protection and comfort after foot surgery or trauma. Adjustable closures expand to accommodate bulky dressings. Square toe design for extra protection of toes, forefoot and bandages. Rocker bottom for a smooth, natural gait. Fits right or left. US Patent D718,024



Sizing	Small	Medium	Large	X-Large
Men's Shoe	10 – 11	11½ – 12½	13 – 14	14½ – 15½
Order #	7500-01MS	7500-01MM	7500-01ML	7500-01MX
Women's Shoe	6 – 7½	8 – 9	9½ – 10½	—
Order #	7500-01WS	7500-01WM	7500-01WL	—

CONTOURED INNER SOLE for Square Toe Post-Op Shoe

7590 Each, SM-XL

Interchangeable with the Powerstep® Square Toe Post-Op Shoe (Model 7500). Sturdy foam with contoured arch support for proper positioning. Provides improved fit, function and comfort of a standard post-op shoe for increased compliance. Deep heel cup cradles the heel for rear foot control. Can be easily trimmed or cored to offload ulcers, incisions or bony prominences. Right or Left specific. US Patent D718,028



Sizing	Men's Small	Men's Medium	Men's Large	Men's X-Large	Women's Small	Women's Medium	Women's Large
Men's Shoe	10 – 11	11½ – 12½	13 – 14	14½ – 15½	—	—	—
Women's Shoe	—	—	—	—	6 – 7½	8 – 9	9½ – 10½
Right Order #	7590-02MS	7590-02MM	7590-02ML	7590-02MX	7590-02WS	7590-02WM	7590-02WL
Left Order #	7590-03MS	7590-03MM	7590-03ML	7590-03MX	7590-03WS	7590-03WM	7590-03WL

PLANTAR FASCIITIS SUPPORT SLEEVE

7012-20 Each, SM-XL

20-30mm Hg compression throughout. Ribbed bands at arch and Achilles tendon offer added support. Ultra-thin, lightweight and breathable toeless design. Measure circumference of foot around arch.



Sizing	Small	Medium	Large	X-Large
Men's Shoe	3 – 5½	6 – 9½	10 – 13	13½ & Over
Women's Shoe	4 – 6½	7 – 10½	11 & Over	—
Arch	5 – 8"	7 – 10"	9 – 13"	11 – 15"
Order #	7012-20SM	7012-20MD	7012-20LG	7012-20XL

DORSAL NIGHT SPLINT

7003-00 Each, REG/LG

Speeds recovery from Plantar Fasciitis and Achilles Tendonitis. Gently stretches the plantar fascia to reduce inflammation and associated pain. Less bulky than traditional night splints. Soft padding throughout and adjustable straps for a comfortable fit. Heel opening prevents migration during sleep. For mild to severe or chronic cases of Plantar Fasciitis. Fits Right or Left. HCPCS L4397 Approved



Sizing	Regular	Large
Men's Shoe	5 – 9	9 – 14
Women's Shoe	6 – 10	10 – 15
Order #	7003-00R	7003-00L

PLANTAR FASCIITIS ARCH WRAP

7100-10 Each, UNIV

Provides targeted support to help relieve arch or heel pain and pain associated with Plantar Fasciitis. Lifts and supports the arch to ease the strain on the plantar fascia. Built-in, contoured gel pad with adjustable strap. Soft materials for longer wear. Fits right or left foot.



Sizing	Universal
Order #	7100-10UN



ULTRASTRETCH® NIGHT SOCK

7001-00 Each, REG/LG

Comfortably stretches the calf and soft tissues of the foot during rest or sleep. Ultra-soft, lightweight and breathable materials. Adjustable straps. Fits Right or Left. Measure around leg at widest part of calf.



Sizing	Regular	Large
Fits Calf	up to 16"	16 – 21"
Order #	7001-00R	7001-00L

ACHILLES SUPPORT

7200-10 Each, UNIV/XL

Cushioning gel pad for targeted compression to support and stabilize the Achilles Tendon. Soft material provides added warmth and compression. Fits right or left foot. Sizing by shoe size.

Sizing	Universal	X-Large
Men's Shoe	8 - 10	10 ½+
Women's Shoe	5 - 9	9 ½+
Order #	7200-10UN	7200-10XL



PATELLA STRAP

7005-10 Each, UNIV

Cushioning gel pad for targeted compression to help offload the strain on the tendon to alleviate pain from overuse injuries. Soft materials with two-way, fully adjustable straps. Universal design fits knees 11-18", measured around leg just below the kneecap. Fits right or left knee.

Sizing	Universal
Order #	7005-10UN



ULTRAFLEX® FOOT ROCKER

7025-50 Each, UNIV

For stretching the plantar fascia, calf muscles, Achilles Tendon, shins, ankles and hamstring muscles.

Stretches and strengthens muscles, ligaments and tendons of the foot and leg. Increases flexibility and strength to help prevent injuries. Non-skid rocker bottom. Stretching exercise guide included.



Sizing	Universal
Order #	7025-50UN

IT KNEE BAND

7008-10 Each, UNIV/XL

Cushioning gel pad for targeted compression over the IT Band to help alleviate pain. Soft materials with fully adjustable closure. Fits right or left knee. Measure around leg just above the kneecap.

Sizing	Universal	X-Large
Fits Knee	Up to 19"	19"+
Order #	7008-10UN	7008-10XL



G2 RECOVERY SOCKS

6120-11
6125-18 Pair, SM-XL

Gradient 20-30mm Hg compression enhances circulation in lower leg for faster recovery. Ribbed bands at Achilles and arch for added support. Knitted with CoolMax® fibers for cooler wear. Seamless toe. Measure around widest part of calf and just above ankle bone.



MEN'S	Small	Medium	Large	X-Large
Men's Shoe	7 & Under	7½ - 10	10½ - 12	12½ & Over
Calf Size	11 - 15"	13 - 17"	14 - 18"	15 - 20"
Ankle Size	7½ - 8½"	8½ - 9½"	9 - 11¼"	11½ - 12½"
BLACK/GRAY	6120-11MS	6120-11MM	6120-11ML	6120-11MX

WOMEN'S	Small	Medium	Large	X-Large
Women's Shoe	4 - 5	5½ - 7½	8 - 10½	10½ & Over
Calf Size	10½ - 14½"	11½ - 15½"	12½ - 17"	13½ - 17½"
Ankle Size	6½ - 8"	8 - 9½"	9½ - 11"	11 - 12½"
BLACK/PINK	6125-18WS	6125-18WM	6125-18WL	6125-18WX



PERFORMANCE CALF SLEEVES

6010-10
6015-10 Pair, SM-XL

Gradient 15-20mm Hg compression for improved circulation while running or during activity. Knitted with CoolMax® fibers for cooler wear. Measure around widest part of calf and just above ankle bone.



MEN'S	Small	Medium	Large	X-Large
Calf Size	11 - 15"	13 - 17"	14 - 18"	15 - 20"
Ankle Size	7½ - 8½"	8½ - 9½"	9 - 11¼"	11½ - 12½"
BLACK	6010-10MS	6010-10MM	6010-10ML	6010-10MX

WOMEN'S	Small	Medium	Large	X-Large
Calf Size	10½ - 14½"	11½ - 15½"	12½ - 17"	13½ - 17½"
Ankle Size	6½ - 8"	8 - 9½"	9½ - 11"	11 - 12½"
BLACK	6015-10WS	6015-10WM	6015-10WL	6015-10WX



With Powerstep®, you are choosing a brand that is recommended by medical professionals nationwide, is proven clinically effective, and has great brand name recognition and consumer satisfaction.

By carrying Powerstep, you can further differentiate your practice with a unique offering not found in retail locations.

Here's how to successfully promote Powerstep® in your practice to maximize revenue:

1 CREATE INTEREST WITH PRODUCT DISPLAYS

Get started with simple organizational solutions created by Powerstep®. We have options to fit your space needs, from compact counter top displays to standing floor spinners and slat wall fixtures.

Powerstep® also offers a professional restocking program to help medical professionals manage and organize office stock of prefabricated orthotics with easy replenishment. Contact your sales representative for details.

2 INTRODUCE PATIENTS TO POWERSTEP

Present patients with the product during the exam and offer Powerstep® as an easy, instant solution to foot pain. Have patients try on styles or use our product style rings to show other suitable options in the Powerstep® offering.



POP Counter Display



Counter top Spinner Display for insoles or footwear



Stock Storage/Display Unit for insoles or footwear



ProTech Style Ring



PowerKids® Sizing Ring

3 **DISPLAY BRANDING PROMINENTLY**

Generate patient interest with our in-office poster in your waiting area or exam room, or with a window cling on the front door or receptionist window. We can customize an ad to your specific practice and/or product focus, as needed.



3' Display Header



Authorized Dealer Window Cling



4 **BE A PATIENT RESOURCE**

Have our informative patient brochures available in waiting areas for patients to learn more about common foot conditions and the benefits of Powerstep®.



5 **PROMOTE YOUR PRACTICE**

Utilize our product images and digital ad templates to promote Powerstep® on your website and social media platforms. We can customize an ad to your specific practice and/or product focus. Also, verify your listing on the www.powerstep.com website zip code practice locator tool to drive business.



Sample social media post



Listing on powerstep.com Practice Locator

To learn more about Powerstep® merchandising and promotional materials, contact your Sales Rep at 888-237-3668.



ABOUT POWERSTEP®

Powerstep® was founded in 1991 by a podiatrist, who was also a foot surgeon and runner. After years of analyzing thousands of casts, it became apparent that there were a few common orthotic shapes that seemed to relieve most patients' foot pain. There was a need for a ready-to-wear orthotic that would effectively relieve foot pain. Over 25 years later, Powerstep® has expanded its offering into new product categories to become a single source for solutions to Plantar Fasciitis and other common conditions.

Powerstep's mission is to provide the most effective products available for preventing and relieving common foot problems. The key to our success has been a combination of effective, innovative products, outstanding customer service and strong brand name recognition by consumers.

Today we are committed to servicing our customers with the quality products they expect from Powerstep®.

SIGN UP* FOR POWERSTEP'S ONLINE CUSTOMER PORTAL

- Easy Online Ordering
- 24/7 Access to Ordering
- View Full Product Catalog
- Access digital assets for merchandising
- Track orders
- Review Past Invoices
- Receive 10% off your first 5 orders placed using the portal*

» **Sign up today at powerstep.com » SIGN IN**

*Available to Powerstep direct dealers only



Customer Service

8930 Global Way
West Chester, OH 45069
Tel: (888) 237-3668
(513) 825-1888
Fax: (513) 825-3381
Email: info@powerstep.com

www.powerstep.com

Check us out on:

